

Labeling Procedure

Each bottle should be labeled as shown on the diagrams below. This should be done before getting the sample, as the writing won't stay on the bottles if they are already wet. The labeling should be done with permanent markers (Sharpies work great!).


Each bottle should contain the following information written legibly in black:

1. Collector ID # - which in our case is always 7488 followed by the:
2. DEP Sequence # - which is different for each site (always 3 digits)
3. Site ID – which is different for each site
4. Date

Small bottles also indicate acid content written in RED.

Each lid indicates collector ID # and DEP sequence #.

Small bottle lids also indicate acid content.



This document was created with Win2PDF available at <http://www.daneprairie.com>.
The unregistered version of Win2PDF is for evaluation or non-commercial use only.